FLEX
Polyclonal Rabbit
Anti-Human
Alpha-1-Fetoprotein
Ready-to-Use
(Dako Autostainer/Autostainer Plus)

Code IS500

ENGLISH	
Intended use	For in vitro diagnostic use. FLEX Polyclonal Rabbit Anti-Human Alpha-1-Fetoprotein, Ready-to-Use, (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. Antibodies to Alpha-1-Fetprotein may be useful for the identification of non-neoplastic and neoplastic liver disease, yolk sac tumors and mixed germ cell tumors (1-3). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Alpha-1-fetoprotein (AFP) is a 70 kDa glycoprotein containing 590 amino acid residues. Cells of the embryonic yolk sac, fetal liver, and fetal intestinal tract synthesize this glycoprotein (4, 5). AFP is a very useful marker for the detection of numerous developmental defects and various pathologies, whether cancerous or not (6). Expression of AFP has been demonstrated in about 44% of hepatocellular carcinomas, and in gonadal and extragonadal germ cells tumors, including yolk sac tumors. AFP has not been detected in normal adult tissue by immunostaining (4). However, traces of AFP are found in normal adult sera and in greater concentrations in maternal and foetal sera and amniotic fluid (7). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
Reagent provided	Ready-to-use polyclonal rabbit antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide.
Immunogen	Alpha-1-fetoprotein isolated from human cord serum.
Specificity	The antibody reacts with human alpha-1-fetoprotein. Traces of contaminating antibodies have been removed by solid-phase absorption with human plasma proteins. In indirect ELISA no significant reaction is observed with normal human plasma. As demonstrated by rocket immunoelectrophoresis, the antibody cross-reacts with the AFP-equivalent protein in cat and wildebeest (8).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation including materials required but not supplied	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). <u>Paraffin-embedded sections:</u> Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is recommended using the 3-in-1 specimen preparation procedure for Dako PT Link. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). Note: After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium. <u>Deparaffinized sections:</u> Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link and following the same procedure as described for paraffin-embedded sections. After staining the slides should be mounted using aqueous or permanent mounting medium. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.
Staining procedure including materials required but not supplied	The recommended visualization system is EnVision™ FLEX, High pH, (Dako Autostainer/Autostainer Plus) (Code K8010). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: AFP (without counterstaining) or AFPH (with counterstaining) The Auxiliary step should be set to “rinse buffer” in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to “none”. This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator’s Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services. Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in “Performance characteristics”. Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin, (Dako Autostainer/Autostainer Plus) (Code K8018). Non-aqueous, permanent mounting medium is recommended.

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- Jacobsen GK, Jacobsen M. Alpha-fetoprotein (AFP) and human chorionic gonadotropin (HCG) in testicular germ cell tumours. A prospective immunohistochemical study. Acta Path Microbiol Immunol Scand Sect A 1983; 91:165-76.
- Kodama T, Kameya T, Hirota T, Shimosato Y, Ohkura H, Mukojima T, et al. Production of alpha-fetoprotein, normal serum proteins, and human chorionic gonadotropin in stomach cancer: Histologic and immunohistochemical analyses of 35 cases. Cancer 1981; 48:1647-55.
- Perrone T, Steeper TA, Dehner LP. Alpha-fetoprotein localization in pure ovarian teratoma. An immunohistochemical study of 12 cases. Am J Clin Pathol 1987; 88:713-17.
- Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 1999. p. 9-10.
- Mizejewski GJ. (-fetoprotein as a biologic response modifier: Relevance to domain and subdomain structure (minireview). Proc Soc Exp Biol Med 1997; 215:333-62.
- Nunez EA. Biological role of alpha-fetoprotein in the endocrinological field: Data and hypotheses. Tumor Biol 1994; 15:63-72.
- Tietz NW. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders Company; 1995. p. 236-38.
- Hau J, Nilsson M, Skovgaard-Jensen H-J, de Souza A, Eriksen E, Wandall LT. Analysis of animal serum proteins using antisera against human analogous proteins. Scand J Lab Anim Sci 1990; 17:3-7.
- Jacobsen GK, Jacobsen M, Henriksen OB. An immunohistochemical study of a series of plasma proteins in the early human conceptus. Oncodev Biol Med 1981; 2:399-410.

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